


AED Series

701 PLUS



AED 701 Plus


Heart Guardian



- Automatic Self-test once a month
- Adult / Child Mode Compatible switch (Energy: 150J ~250J for Adult / 50J for Child)
- 4.3" High Resolution Color LCD display
- Size: 222(W) x 280(D) x 75(H)mm
- Weight (with Battery) : 2.3kg (5lb)


AED 701

Heart Guardian



- Automatic Self-test once a month
- Adult / Child Mode Compatible switch (Energy: 150J ~250J for Adult / 50J for Child)
- Visual Guidance with LED Lights
- Size: 222(W) x 280(D) x 75(H)mm
- Weight (with Battery) : 2.3kg (5lb)

P-701 (701Series Pads)



- Type : Pre-Gelled / Self-Adhesive / Disposable Non-Polarized Defibrillation Pads Non-Sterilization
- Shelf-Life : 2 years from Date of Manufacturing
- Size & Length : 110 x 120mm / 1.8m
- Pouch Size : 170 x 250mm

BT-701 (701Series Battery)



- Type : LiMnO2(12V / DC6.4A) / Non-Rechargeable
- Weight: 285g
- Shelf-Life : 4years for Standby Mode (Can be reduced by temperature and Humidity condition during storage)

B7 (701Series Semi-Rigid Carrying Case)



Closed

Open