



EVOTEC
EvoCare®

Product Catalogue

 **ceasa**
THE CLINICAL ENGINEERING ASSOCIATION
OF SOUTH AFRICA
corporate member

SAHPRA

License No:
00000598 MD

658 Alouette Street
Elarduspark
Pretoria

P.O. box 4333
Rietvalleirand
0181

Head Office: Tel : 012 345 1176 / 4738 / 6421
Fax : 012 345 4381 or 086 550 4126
E-mail: sales@evotec.co.za

Cape Town : 082 447 3175 **George**: 071 426 1012 **Johannesburg** : 060 783 6163
Pretoria : 073 358 7258

www.evotecmedicalsupplies.africa

A	Accessories	1, 2, 5, 8, 14	J	Jaundice detector	11
	Adaptor cables		L	Lamps - ENT , Theatre , Various	1
	- Electrosurgical	1		Laryngoscopes	6
	- IBP/SpO2	2		Latex free cuffs	4
	- ECG	2		Lights - Theatre / Examination	13
	AED Pads	1		Light sources	14
	AED's	5		Lung Function	14
	Anaesthetic	3	M	Mattress	24
	Aneroids	24		Measuring tape	19
	Autoclaves	3		Monitoring Accessories	1, 2
	AVALON Transducer repairs	8		Monitoring straps - CTG	8
B	Batteries	1		Monitors - Vital Signs	16
	Baby cart	11		Monitors - Handheld / Finger	16
	Baby Doppler	11		Monitors - ICU & Theatre	16, 17
	Baby scales	11, 19		Monitors - Fetal	15
	Bags - feeding	12		Monitors - oxygen	11
	Bathroom scales	19	N	Nebulisers	18
	Blenders - Oxygen	18		Neonatal warmers	10
	Blood pressure Aneroids / Arm / Wall	24		Neonatal Phototherapy	11
	Blood pressure cables	2		NIBP Cuffs / hoses / connectors / bulbs	1
	Blood warmers	4	O	Ophthalmoscopes	6
	Blood pressure meters	4		Otoscopies	6
	Bodoc seals	9		Oxygen Blenders	18
	Bottles - suction	2		Oxygen Cells	18
	Bottles - Humidifier	9		Oxygen Cylinders / bags / trolleys	18
	Bullnose regulators	8		Oxygen generators	18
C	Cables - Patient / recorders	1, 2		Oxygen Hospital Plant generation	18
	Cardiac monitors	16, 17		Oxygen monitors	11
	Cells - Oxygen	18	P	Pacemaker module	17
	Cleaners - Ultrasonic	3		Paper	1, 24
	Connectors NiBp / ECG	1		Patient cables	1, 2
	CTG - Gel	8		Patient transfers	24
	CTG monitors	15		Peak flow meters	14
	Cuffs - Blood pressure	1, 4		Phototherapy	11
D	Defibrillators & PADS	5		Plaster saws	24
	Diagnostic sets	5, 6	R	Rebreathing bags	18
	Diagnostic wall units	5		Recorders	7
	Disposables	25		Regulators	8, 9
	Diverters - Peep valve	18		Resuscitators	18
	Dopplers - Baby / Fetal / Vascular	11, 23	S	Seals & O-rings	9
	Drivers syringe	12		Service & repairs	25
E	ECG Recorders	7		Scales	19
	Electrosurgical	7		Silicone tubing	2
	Electrical safety analyser	21		SpO2 probes	2
	Examination Lights	13		Sterilizers	3
F	Feeding pumps & Bags	12		Stethoscopes	24
	Fetal Dopplers	23		Suction filters & bottles	2
	Fetal transducers	8		Suction machines	20
	Fetal Monitors	15		Syringe drivers	12
	Fetal transducer repairs	8	T	Temperature probes / sensors / Infrared	2, 24
	First aid kits	24		Test equipment	21
	Flow meters	8, 9		Transcutaneous Jaundice Detector	11
	Fluid warmer	4		Transducers - Fetal -New and refurbished	8
	Furniture	24	U	Ultrasounds / Printers / Paper	22
G	Gauges - Regulators	9		Ultrasonic Cleaners	24
	Gel - CTG	8		UV Germ inspection box and gel	3
	Gel - Ultrasound	22	V	Vaporizers	3
	Generators - Oxygen	18		Various Items	24
	Germ UV inspection box , gel and powder	3		Vascular dopplers	23
	Glucometers	24		Vital signs monitors	16
H	Handles - Diagnostic	6	W	Warming Cabinets	4
	Humidifier - Anaesthetics	3		Welch Allyn	6
	Humidifier - Flow meters	9	X	X-Ray Viewing Boxes	14
I	Incubator - Accessories	1			
	Incubators	9, 10			
	Infant resuscitators	18			
	Infant warmers	10, 11			
	Infusion pumps	12			
	Infusion pump docking station	12			
	Inspection box - Germ UV	3			



AED Pads

With or without leads

Cardioline	CU Medical
Defibtech	Marquette
Hewlett Packard (HP)	
Instramed	Mediana
Medtronic Physiocontrol	
Metsis	Nihon Kohden
Philips	Schiller
Welch Allyn	Zoll



Batteries

Button	Lead Acid
Lithium	Custom Packed

Philips , Nihon Kohden , GHealth
BLT, Sinohero, 3F , Adecon, GE,
Drager, Medtronic , etc



Electrosurgical

Pencils—Reusable &
disposable Pads— Reusable &
disposable Adaptor cables
Pens & LED pens
Forceps
Suction accessories



Incubator and Warmers

Filters
Sleeves and cuffs
Hydrogel Probe covers
Mattresses
Temperature probes



Lamps

Heine , Welch Allyn , Riester ,
Kawe, Hanualux, Ge Lighting,
Osram, Dr Fisher, Zumax

Infrared Ultra violet Theatre
ENT, Examination , etc



NiBp Cuffs

Super tough and super soft cuffs

Neo Natal
Child
Adult
Large Adult
Thigh
1 and 2 tube

A-LOCK
E-LOCK
SUPERSOFT



NiBp Connectors

Luer	IRIS
Barbed	Philips
Marquette	Siemens
Snap Luer	Redel
Dinamap	Monitor
Nihon Kohden	Cuff



NiBp Hoses

Luer	IRIS
Barbed	Philips
Marquette	Siemens
Snap Luer	Redel
Dinamap	
Nihon Kohden	



NiBp Various


Bulbs
Valves
Coiled tubing 4 and 8 foot
Hoses (without connectors)



Patient Cables

3 lead - 1 piece & 2 piece
5 Lead - 1 piece & 2 piece
10 Lead - 1 piece & 2 piece
Multi Cables

All Monitor Manufacturers



Paper

Recorders
CTG
Ultrasound
Monitors
Defibrillators

Pictures are for illustrative purposes only and may vary from actual product.



Patient Cable Adaptors

Multi electrode adaptors
Banana to Clip
Banana to clip with lead



SpO2 Probes - 3m

For Adults , Neonates & Children
3meter complete from Monitor to patient
Clip , Silicone Soft , Wrap, Ear
ALL MANUFACTURERS



SpO2 Probes - 0.9m

For Adults, Neonates, Children
Needs Adaptor cable
Clip , Silicone Soft, Wrap, Ear
ALL MANUFACTURERS



SpO2 / IBP Adaptor Cables

Adaptor Cables for connection to 0.9m SpO2 probes
Adaptor cables for invasive blood pressure transducers
ALL MANUFACTURERS



SpO2 - Various

Connectors
Cable
Clips , boots, silicone soft sleeves, ear Clips

**WE DO REPAIRS ON SPO2 .
(Certified test equipment)**



Silicone Tubing - Biomedical Grade

100 m rolls to 1 meter lengths
Consistent inner / outer diameters
High grade at distributor prices



Suction Filters

Replacement Suction Bacterial filter for all suction units



Suction Bottles

1L - Bottles
2L - Bottles
4L - Bottles

With or without liners
Perspex or glass



Temperature Probes

Reusable and disposable
YSI 400 , YSI 700 , 10K series

Siemens , Philips, Mindray ETC

Monitoring and Incubator sensors
FOR ALL MANUFACTURERS



**Various accessories not listed
PLEASE CALL**

Pictures are for illustrative purposes only and may vary from actual product.



Anaesthetic Units

Morpheus E units from Europe
Highest standard
3 gas mixing
2 Vaporizers
Electronic ventilator
12" TFT screen



Anaesthetic Units

Morpheus LT units from Europe
Highest standard
3 gas mixing
2 vaporizers
Electronic ventilator with
12" TFT colour screen



Siare Vaporizers

Sevoflurane
Fluo- Halothane
Ethrane
IsoFlurane
DesFlurane



Anaesthetic Options

Spirometry gas monitor
Patient circuit support arms
Scavenger for Anaesthetic gas
Tracheal Venturi aspirator
Etc

Autoclaves & Sterilizers



Autoclaves 16 / 24 L

CLASS B MA series
Pre - Vacuum
Dry - Vacuum
Microprocessor and Temperature
control
Special Liquid sterilization program




Autoclaves 16 / 24L

CLASS S FLASH
Unwrapped / Wrapped
Fully Automatic
Pre heat an Dry functions
Auto Water Add function



Autoclave Vertical

50 to 300 L capacity
Fully automatic
Unwrapped / Wrapped
Auto add water function
2KW heater power



Sterilizers & Cleaners

UV sterilizers 8 / 15 / 30W



Ultrasonic cleaners / baths



UV Germ inspection box

Used together with UV germ
lotion and UV germ powder to
train aseptic techniques.
Comes with optional carry case



GLO GERM Gel & Powder

Glows under UV light
For training of aseptic techniques
Handwashing, Gloving, Gowning,
Masking, Isolation, Surface
cleaning
To be used with UV box or UV
torch

Pictures are for illustrative purposes only and may vary from actual product.



Blood / Fluid warmer

Constant temperature 37 ° C
Auto Temperature regulated
Small and portable
On / Off
Ease of use

Evo-bw



Blood / Fluid warmer

Microprocessor controlled
28 to 42 °C
Audio and Visual alarms
Continuous operation
Waterproof

BF1000



Blood warmer Single

Microprocessor controlled
TFT LCD display
Liquid sense detector
6 warming modes
Alarms for hi/low, no liquid, fault sensing
9L / H

BFW-1020S



Blood warmer Dual

Microprocessor controlled
TFT LCD display
Liquid sense detector
6 warming modes
Alarms for hi/low, no liquid, fault sensing
18L / H

BFW-1020D



Warming Cabinets

Microprocessor controlled
TFT LCD display touch screen
10 L, 30L, 50L

Blood pressure meters



BP Meter Handheld

Complete with aneroid
Bulb with control valve
Pouch

Latex free cuffs available



BP Meter Aneroid Wall

Complete with coiled tubing
Bulb with control valve
Swivel
Cuff

Latex free cuffs optional available



BP Meter Aneroid Mobile

Complete with coiled tubing
Bulb with control valve
Mobile trolley
Cuff

Latex free cuffs available



BP Meter electronic Desk

Fully automatic measurement of
Systolic, diastolic and pulse
Uses FUZZY logic
80 memories and average test results
Automatic power off



BP Meter electronic Arm

Automatic measurement of Sys-
tolic, Diastolic and pulse
Low noise pump
Auto power off
History memory



LATEX FREE CUFFS Available

A-Lock E-Lock Super Soft

Neonatal Child
Adult Large Adult
Thigh


1 and 2 tube

Pictures are for illustrative purposes only and may vary from actual product.



AED HR501

1 2 3 defibrillation
Pacemaker detection
Motion detection
Lightweight
Switch selection between
Adults and Paeds - no pad
swop



AED HD701 with screen

Biphasic
LCD screen
2 - 270 Joule
Switch selection for Adults/
peads - No need to change
pads



Defibrillators

Biphasic
AED
SDP Sudden death prevention
12 lead ECG
NiBp , SpO2, EtCo2 , Pacing
Recorder
User replaceable Battery



AED PADS

Adult and Paediatric pads
Standard and Lead Out pads
Radiotranslucent
Paediatric reduced energy

For ALL manufacturers

Diagnostic

WelchAllyn

zumax

We sell Welch Allyn and
Zumax Diagnostic equipment



Diagnostic Set Budget

Black with 2 spare bulbs



Diagnostic Set ZUMAX

Highest quality
Precision optics and
vacuum lamps



Wall Diagnostic Set ZUMAX / Welch Allyn

Wall mounted modular ENT
High precision optics
Memory of light settings



485 Wall mount trans- former ZUMAX / Welch Allyn

Wall mount transformer with
extension unit and power
cable



Welch Allyn Wall Diagnostic Connex

Wall mounted ENT set
High precision



Welch Allyn Diagnostic sets Various

Pictures are for illustrative purposes only and may vary from actual product.



Ophthalmoscope ZUMAX

3.5V Coaxial
3.5V PRO-Coaxial
3.5V LED

With and Without handles



OTOSCOPE ZUMAX

3.5V LED
3.5 Fiber

With / without handle (K-Handle or Rechargeable)



Specula ZUMAX

Sizes 2, 3, 4, 5, 9



Laryngoscopes

Conventional & Fibre Optic
2.5V & 3.5V LED



Desktop Charging unit ZUMAX

DT801—Single charger
DT802—Double charger



Desktop Charger for K-Handle ZUMAX

DT407 - Single
DT408 - Double



Handles ZUMAX

490-0010 3.5V Metal Li-Ion re-chargeable
490-0015 3.5V Metal with USB rechargeable
490-0019 Metal for C size batteries rechargeable



Diagnostic sets ZUMAX

Rechargeable diagnostic sets

Ophthalmoscope
Otoscope or
Retinoscope



Diagnostic Sets ZUMAX

Different sets to your specifications



WELCH ALLYN

EVOTEC also supplies Welch Allyn products with full after sales service

Pictures are for illustrative purposes only and may vary from actual product.

Recorder - 1 and 3 channel Interpretive available

1, 2, 3 c (- 2) channel printout
3.8inch LCD display with waveform display for 3 / 6
12 leads
Auto measurement and 128 exam memory
USB and data interface



Recorder - Holter

PM85 for 24hour dynamic monitoring
10 template user-defined modules
Flexible analysis channel
Flexible Atrial fibrillation Analysis
Pacemaker analysis
Case management
Sleep apnea analysis

Recorder - 6 channel Interpretive

3, 3+1, 3+6, 6 , 6+1 printout
35 inch LCD display with wave-
form display
Simultaneous 12 lead acqui-
sition.
USB to PC



Recorder - 12 channel Interpretive

6/ 6+1 / 12 channel printout 6
inch movable screen 12 lead
display , Touch screen
Simultaneous 12 lead acqui-
sition. Analog to digital
conversion at patient

Recorder - 12 channel ECG1200G Telemedicine

6/ 6+1 / 12 channel printout High
Res LCD screen Simultaneous
12 lead acqui-sition. Waveform
display Touch screen, wifi, 3G,
cloud Management system ,
Tools etc



Recorder - 12 channel Interpretive ECG1200G

6/ 6+1 / 12 channel printout High
Res LCD screen Simultaneous
12 lead acqui-sition. Waveform
display
USB to PC. Tools Menu .
Management system

Electrosurgical

OBS-50 High Frequency Coagulator

Widely used in plastic surgery,
dermatology, ENT,
ophthalmology, gynecology
proctology and urology.



Electrosurgical CCS-2001

0 - 300 W
Single and dual frequency
Cut , Coag , bipolar
Ground plate detection

Electrosurgical HV-300 plus and HV-400 plus

0 - 300 W and 0-400
Nine working modes Monopolar:
Pure cut, blend1, blend2, blend 3,
spray coag, forced coag and soft
coag. Bipolar: Bipolar coag,
bipolar cut REM . Endo-safe and
Seal-safe



Electrosurgical CCS-2000AI LCD

0 - 400 W
9 output modes
Bipolar and monopolar
Under Water
REM warning system
Remote Control

Pictures are for illustrative purposes only and may vary from actual product.



Transducer - TOCO EX STOCK

Bistos Bionics
Corometrics Contec
MARS Marquette
Edan Philips
GE—Medical HP
Unicare others not listed



Transducer - US EX STOCK

Bistos Bionics
Corometrics Contec
MARS Marquette
Edan Philips
GE—Medical HP
Others not listed



Transducer repairs

Philips AVALON
Transducer repairs !!!!
Refurbish !!!

**New casings, cables,
strain reliefs, connectors**

Other manuf. Also available



CTG Gel & Fetal Monitor straps

Conductive gel for use with ultra-
sounds and CTG units

Fetal monitoring straps for ulti-
mate comfort— both TOCO and
US

Flow-meters and regulators



Flowmeter - Single Tube type

0 - 15 L / min
With "Heyer" connector
Flow outlet included

Silver or white is for Oxygen
Black = Air Green = Carbon
Dioxide



Flowmeter - Double Tube type

0 - 15L / min
With Heyer connector
Flow outlet included

Silver or white is for Oxygen
Black = Air Green = Carbon
Dioxide



Flowmeter - Single Click type

0 - 15 L / min
With "Heyer" connector
Flow outlet included

Silver or white is for Oxygen
Black = Air Green = Carbon



Flowmeter - Double Click type

0 - 15 L / min
With "Heyer" connector
Flow outlet included

Silver or white is for Oxygen
Black = Air Green = Carbon



Regulator - Bullnose

0 - 15L/min (others
available) High pressure

With convenient lock wheel (no
spanners needed)
Bullnose to Pin-Index adaptors
available



Regulator - Pin Index

0 - 15L/min (others available)
High pressure

Pictures are for illustrative purposes only and may vary from actual product.



Flow Output

Robust
High impact plastic outlet for flowmeters and regulators



Humidifier Bottle

Can be direct connected to the flowmeters
Robust and lightweight



Regulator Seals and O-rings

Special seals for Pin Index regulators
O-rings for Bullnose regulators



Regulator Gauges

Spares - regulator gauges for pin-index and bullnose click type regulators

Incubators



Incubator BB-100S

Microprocessor Based Servo controlled
Air mode temperature system
Large LED indication
Five Alarm functions Visual and audible
Bed Inclination adjustable
Humidity adjustable in 2 steps



Incubator BB-100T

Microprocessor Based Servo controlled
Air mode temperature system
Large LED indication
Five Alarm functions Visual and audible
Bed Inclination adjustable
Humidity adjustable continuously
MORE ...




Incubator BB-200S

Microprocessor Based Servo controlled
Air and skin mode temperature system
Air , Skin temp, Time and heating power with Large LED indication
Five Alarm functions Visual and audible
Bed Inclination adjustable
Humidity adjustable continuously
MORE




Incubator BB-200T

Microprocessor Based Servo controlled
Air and skin mode temperature System
Humidity display system
Air , Skin temp, Time and heating power with Large LED indication
Five Alarm functions Visual and audible
Bed Inclination adjustable
Humidity adjustable continuously



Incubator BB-200L

Microprocessor Based Servo controlled
Air and skin mode temperature System
Humidity display control system
Air , Skin temp, Time and heating power with Large LED indication
Five Alarm functions Visual and audible
Bed Inclination adjustable MORE ...



Incubator BB-300S

Microprocessor Based Servo controlled
Air and skin mode temperature System .
Humidity and Oxygen control system
Air , Skin temp, Time and heating power with Large LED indication .
Five Alarm functions Visual and audible .
Bed Inclination adjustable
Humidity adjustable continuously .
Double hood insulation. MORE

Pictures are for illustrative purposes only and may vary from actual product.



Incubator BB-300T

Microprocessor Based Servo Controlled. Air and skin mode temperature system . Humidity and oxygen control Air , Skin temp, Time and heating power with Large LED indication . Five Alarm functions Visual and audible . Bed Inclination adjustable . Electric height adjust . Double hood Humidity adjustable continuously .



Incubator BB-300 L

Microprocessor Based Servo Controlled. Air and skin mode temperature system . Humidity and Oxygen control system Weighing system . Air , Skin temp, Time and heating power with touch screen and graphs. Five Alarm functions Visual and audible . Bed Inclination adjustable . Electric height adjust . Humidity adjustable continuously . Double hood insulation. MORE



Incubator Transport BT-100

Microprocessor servo controlled Air and Skin mode temperature AC and DC power—DC 12V or 24V SET, AIR and BABY temperature displayed. Battery indication Double wall hood with Side door to pull out bed. Height adjustable . Oxygen cylinder and oxygen supply system . Adjustable observation lamp



Options for Incubators and warmers

- Tray
- Electronic height adjuster
- Double hood
- Baby scale
- X-ray cassette

Neonatal Warmers



BN-100 Standard

Microprocessor Controlled servo Pre-Warm , Manual or automatic Audible and Visual alarms Apgar Timer. X-Ray cassette tray Heater Overhead can be swivelled 90 degrees. Incline of bed is adjustable. Adjustable LED observe light . Transparent protector



BN-100 Top Grade

Microprocessor servo controlled Prewarm, manual and automatic Self test, Apgar timer, X-ray cassette, Bed incline adjustable, Overhead crystal heater adjustable. LED observe light adjustable. Built in low pressure aspirator. LED phototherapy. Many options available



BN-100 Luxurious

Microprocessor servo controlled Prewarm, manual and automatic Self test, Apgar timer, X-ray cassette, Bed incline adjustable, New design carbon steel heater. LED observe light adjustable. Large bed. Built in low pressure aspirator. LED phototherapy. Oxygen supply. Many options available



BN-200 Luxurious

Microprocessor servo controlled Prewarm, manual and automatic Set, Skin and peripheral temperature 8" colour touch screen. Integral timer Self test, Skin temp trend. X-ray cassette, Bed incline adjustable, Electric height adjustment. New design carbon steel heater. LED observe light adjustable. Large bed. Built in low pressure aspirator. LED phototherapy. Oxygen supply.



BN-100A Standard

Microprocessor controlled Audible and Visual alarms Transparent protector Second thermal switch off Skin sensor drop-off sensor to avoid over temperature



BN-100A Top Grade


Microprocessor Controlled Pre-Warm , Manual or automatic settings . Audible and Visual alarms Apgar Timer. X-Ray cassette tray Heater Overhead can be swivelled 90 degrees. Incline of bed is adjustable. Adjustable observe light . Transparent protector Skin sensor drop off protection

Pictures are for illustrative purposes only and may vary from actual product.



BN-200

Microprocessor based Servo Controlled system. Manual, skin and mattress settings . Audible and Visual alarms. Apgar Timer. X-Ray cassette tray. Heater Overhead can be swivelled 90 degrees. Incline of bed is adjustable. Adjustable LED observe light .
Integrated resuscitation unit
Integrated Suction
Oxygen supply system & air temp



BL-30D/ BL100

LED phototherapy
Button control
High intensity
No Noise
Integral therapy timer
Count down timer
30D -Top mount of incubator
100—under bed—auto radiance adjustment



Phototherapy BN-300

3.5" tFT touch control screen
Light intensity adjustable 0-100%
32 pcs Super LEDs and 8 pcs white examination LEDs
Real time clock and date
0-99 min treatment timer
Can be disconnected from trolley for separate use.



Phototherapy BL-50 LED

Button controlled LED radiating elements
High Intensity
Equal light distribution
Adjustable head and height
Long lifespan
Integral therapy timer

Neonatal-VariouS



Baby Cart

Adjustable height
Trendelenburg adjustment
Available in Green , Pink and Blue
High quality with aluminium column



Baby Respiratory Monitor

Alarm if no movement detected for 20 sec
Alarm if breath slows down to less than 10 breaths / min
Compact and portable—no mains cables



Baby Doppler Sonoline C

Fetal Heart rate, bar graph and waveform
Colour screen
Built in speaker , headphone socket . MORE



Oxygen Monitor MAX-02

Microprocessor Controlled
Continuous monitoring of Oxygen
Large display
Flexible alarms



Transcutaneous Jaundice Detector

Dynamic clinical examination
Of neonate jaundice
LCD display



Baby Scales

Various baby scales - electronic
Please call for specifications

Pictures are for illustrative purposes only and may vary from actual product.



ZNB-XD

Curvilinear peristaltic
Uses any IV set
Flow rate 1- 1100ml/h
Purge and Bolus
Volume 1-9999 ml
Infusion modes : ml/h and drop/min
KVO , Alarms, Battery backup



ZNB-XK

Large display
Curvilinear peristaltic
Uses any IV set
Flow rate 1- 1300ml/h
Purge and Bolus
Volume 1-9999 ml
Infusion modes : ml/h and drop/min
KVO , Alarms, Battery backup



ZNB-XA

Curvilinear peristaltic
Uses any IV set
Flow rate 1- 1100ml/h
Purge and Bolus
Volume 1-9999 ml
Infusion modes : ml/h and drop/min
KVO , Alarms, Battery backup
Occlusion adjustable



ZNB-XAII

Large display
Curvilinear peristaltic
Uses any IV set
Flow rate 1- 1300ml/h
Purge and Bolus
Volume 1-9999 ml
Infusion modes : ml/h and drop/min
KVO , Alarms, Battery backup

Infusion lines for ALARIS . Universal lines also available

Syringe Drivers & Feeding



KL-602

Syringe size: 10, 20 , 30 , 50ml
Built in brands: BD, B Braun ,
Terumo, 2 brands user defined
Volume : 0.1-9999ml
Bolus , KVO, Syringe ID, Occlu-
sion, mute, system memory
Alarms



KL-702

Syringe size: 10,20,30,50/60ml
Built in brands: BD, B Braun ,
Terumo, 2 brands user defined
Volume : 0.1-9999ml
Bolus , KVO, Syringe ID, Occlu-
sion, mute, system memory
Alarms




Feeding Pump EA113

Double mode: Thick and normal
Dedicated feeding disposables for better accuracy Built in double microprocessor for best
reliability
Visual audible alarms
Self priming pump door. Internal lights and transparent door for easy identification
Bubble & occlusion detection with embedded sensors
Data confirm key for better accuracy



Feeding Bags

For Feeding pump EA113
Evotec also supplies various oth-
er feeding bags



Docking Station

More pumps, one system
4-in-one syringe pump + 2 infusion
pumps
1 unified power cord

Pictures are for illustrative purposes only and may vary from actual product.



LED Lights Shadowless

Cluster 5+5 160 000 Lux
Cluster 5+3 160 000 + 100 000 LUX
Cluster 3+3 100 000 + 100 000 LUX
Cluster 5 160 000 LUX
Cluster 3 100 000 LUX wall & ceiling



LED Lights Shadowless with colour adjustment

Cluster 5+5 160 000 Lux
Cluster 5+3 160 000 + 100 000 LUX
Cluster 3+3 100 000 + 100 000 LUX
Cluster 5 160 000 LUX
Cluster 3 100 000 LUX wall & ceiling



LED-0705 series

Advanced colour temperature Control
LCD display control system
German Light source
Precise illumination
Compact
Adjustable to 180000 lux



YD02-LED3s

Emergency mobile—shadowless with colour adjustment
140 000 lux
Colour reduction
57 LED lamps
Light field : 150 to 200 mm



LED-3E

Emergency Mobile - shadowless with colour adjust
100 000 LUX
Emergency time 3 hours
Spot diameter 200-350 mm



200 S / E - LED

LED Cold light operating lamp
25 000 LUX
220V A with battery backup
Spot diameter 150 mm

Lights - Examination



YD-01-1SA LED

Examination Light—Intensity adjust.
20 000 LUX
120mm spot diameter
LED— Super bright lumitron



YD-01-1SE LED (Backup)

Examination Light—Intensity adjust.
20 000 LUX
120mm spot diameter
LED— Super bright lumitron
6 hours backup time



YD01A (LED)

Examination light - intensity adjust
Mobile - Gooseneck
30 000 LUX
Spot diameter 100mm
LED - Superbright Lumitron



YD01W (LED)

Examination light - intensity adjust
Wall - Gooseneck
30 000 LUX
Spot diameter 100mm
LED - Superbright Lumitron

Pictures are for illustrative purposes only and may vary from actual product.

Lights - Examination



QS - Q3 LED

Q3 = mobile Q3W = Wallmount
15 000 LUX
Colour temperature 6000 K
LED power 3W
Voltage 220V AC 50Hz



QS - Q5 LED

Q5 = mobile Q5W = Wallmount
Luminous Flux
Colour temperature 6000 K
LED power 3W
Voltage 220V AC 50Hz
Focus control



1082 / 2082 / 3082 LED

1082- Mobile 2082-Wall 3082 Desk
12 000 LUX
Colour temperature 6000K
LED Power 4W
Voltage 220V 50 Hz



Light Source LED / Halogen

Various Lightsources available
Halogen
LED



X-Ray boxes

Single , Double, Triple and quad
boxes
Standard and light regulated
typetube type



LED - Headlight

LED 3W
320000 LUX
DC / AC power supply
4 hours battery time

Lung Function



SP10W Spirometer

FVC, FEV1, FEV1%, PEF, FEF25
FEF75, FEF 2575
Measured and predicted value
Wireless transmit
Flow rate - Volume chart
Volume - Time chart , Trend Chart
Data memory , delete, upload and
review
Rechargeable Lithium battery




SPM-W Spirometer

FVC, VC, MVV + 32 parameters
Measured and predicted value
Respiration waveform in realtime
Flow-volume loop , volume-time
chart, Trend chart
USB or Bluetooth data transmit
Rechargeable lithium battery
Power and measure indicator



Peak Flow Meter

LCD Display with low power con-
sumption
Measures PEF, FEV1, FVC
Health status indication
Calibration function
USB upload
PC software with trend chart view,
range display, time sync and pdf
report output
Low power indicator



Mouth Pieces & noseclips

Disposable mouthpieces
Noseclips

Pictures are for illustrative purposes only and may vary from actual product.

MARS - B



Feature	Fetal Parameter					Maternal Parameter	Printer Paper Width		Auxiliary Device			Advanced Fischer NST Report
	Single FHR	Dual FHR	T O C O	AFM/ MFM	Signal Overlap Alarm	ECG, NIBP, SpO2, PR, TEMP, RESP	112 mm	150 mm	Lithium battery	Touch Screen/ Keyboard / Mouse/ Central Monitoring Station/ HDMI	Fetal Stimulator/ Waterproof Probe/ Wall-mount/ Trolley	
Config.A	✓	✗	✓	✓	✗	✗	✓	Opt	✓	Opt	Opt	Opt
Config.B	✓	✓	✓	✓	✓	✗	Opt	✓	✓	Opt	Opt	Opt
Config.C	✓	✗	✓	✓	✓	✓	✓	Opt	✓	Opt	Opt	Opt
Config.D	✓	✓	✓	✓	✓	✓	Opt	✓	✓	Opt	Opt	Opt

MARS - K



Feature	Fetal Parameter					Maternal Parameter	Printer Paper Width		Auxiliary Device			Advanced Fischer NST Report
	Single FHR	Dual FHR	T O C O	AFM/ MFM	Signal Overlap Alarm	NIBP, SpO2, PR	112 mm	Lithium battery	Touch Screen/ Keyboard / Mouse/ Central Monitoring Station	Fetal Stimulator/ Waterproof Probe/ Wall-mount/ Trolley		
Config.A	✓	✗	✓	✓	✗	✗	✓	✓	Opt	Opt	Opt	
Config.B	✓	✓	✓	✓	✓	✗	✓	✓	Opt	Opt	Opt	
Config.C	✓	✗	✓	✓	✓	✓	✓	✓	Opt	Opt	Opt	
Config.D	✓	✓	✓	✓	✓	✓	✓	✓	Opt	Opt	Opt	

T10 / T12



10.2" LED display / 12.1" LED display
 Flip screen, multi angle, backlight
 Wide beam multi crystal waterproof probe. Stable signal detection
 Flat surface waterproof TOCO transducer
 60 hours storage and replay memory.
 Fisher analyses results
 High Capacity lithium battery
 Options : Twins , Maternal SpO2 / NiBp / Temp.
 Wireless and wired networking

Fetal Monitor accessories

Recording paper



Ultrasound Gel



CTG belts



Trolleys



Pictures are for illustrative purposes only and may vary from actual product.



H5 Vital signs

7 inch TFT colour screen
H5 - ECG, SpO2, Nibp, Temp, Resp
H5A - SpO2, PR
H5B - SpO2 , ETCO2, PR, Resp
H5C - SpO2 , Nibp, PR
H5D - ECG, Temp, Resp, PR
H5E - SpO2, NiBp , Temp, PR
H5F - ETCO2, Resp
H5G - SpO2, NiBp, ETCO2, Resp,PR
H5H - NiBp, SpO2, ECG, RESP, PR



H6

8.4" inch TFT Colour Screen (large display of parameters)
Graphic and Tabular Trends
Parameters : ECG, Respiration, NiBp, Temperature and SpO2 as standard
Adult , Paediatric , Neonatal settings
Battery backup : li-ion
Internal power supply
Options : Touch screen , ETCO2 , Massimo spo2 , IBP etc.



H8

12 inch TFT Colour Screen (large display of parameters)
Graphic and Tabular Trends
Parameters : ECG, Respiration, NiBp, Temperature and SpO2, OxyCRG as standard
Adult , Paediatric , Neonatal settings
Battery backup : li-ion
Internal power supply
Options : Touch screen , Co , Massimo SpO2 , IBP , gasses etc.



H9

15 inch TFT Colour Screen (large display of parameters)
Graphic and Tabular Trends
Parameters : ECG, Respiration, NiBp, Temperature and SpO2 , oxyCRG as standard
Adult , Paediatric , Neonatal settings
Battery backup : li-ion
Internal power supply
Options : Touch screen , Co , Massimo SpO2 , IBP , gasses etc.



H100 handheld

2.8 inch TFT Colour screen
Adult, Pediatrics and neonates
H100C - ETCO2 , PR, Resp
H100S - SpO2 , PR
H100SC - ETCO2, SpO2, PR, Resp
Alarms—visual and audible



OPTIONS

Touch Screen
Invasive BP—2 channels
CO2
Multi Anaesthetic gasses (IRMA) Swiss
IOC module - Index of Cerebral Central Station
Massimo SpO2
Recorder
Brackets or trolleys



Vital Signs PC900 PRO

Portable vital signs monitor
PC900PRO SN - SpO2 + NiBp
PC900PRO SNT - SpO2 + NiBp + Temp
PC900PRO SNET—SpO2 + NiBp + ECG + Temp
Backup battery with internal power supply



Vital Signs MD2000

MD2000 A - SpO2
MD 2000B - SpO2 + NiBp
Backup battery with internal power supply



Vital Signs Handheld

Portable vital signs monitor
PC303— NiBp + SpO2 + Temp
PC304—NiBp + SpO2 + Temp + GLU
PC305—NiBp + SpO2 + Temp + GLU + ECG
All unit incorporate Bluetooth
Internal rechargeable batteries



SpO2 Finger Clip

Display : SpO2
SpO2 waveform
Heart Rate - HR
Battery strength
Colour display
Alarms for SpO2 and HR
Size: 60mm x 35mm x 32mm
Weight : 50 grams

Pictures are for illustrative purposes only and may vary from actual product.

HT – 6 MODULAR



Features

- 12.1" colour TFT display
- Support touch screen
- Auto ID of plug-in modules
- Multi display
- Standard: ECG, NiBp, SpO2, Temp, Resp
- Optional: ETCO2, IBP, CO, EEG, Anaesthetic gasses, Depth of Anaesthesia, Mini Host, Pacemaker

Data & Storage

- 72 hour trend , optional thermal recorder, real time ecg and spo2 waveform
- ST analysis—13 sets , Drug dose calculation—15 sets
- VGA output, BNC connection , USB, SD card
- Wired / wireless Central monitoring system

Mini Host & Pacing

- HT10 MINI HOST can transfer data from different units
- ST analysis—13 sets , Drug dose calculation—15 sets
- VGA output, BNC connection , USB, SD card
- Wired / wireless Central monitoring system

HT – 8 MODULAR



Features

- 15" Colour TFT display
- 1024*768 resolution
- See HT 6

HT – 9 MODULAR



Features

- 17" Colour TFT display
- 1280*1024 resolution
- See HT 6



HT-10 Mini Host

- Intelligent module to be used separate or with HT line of monitors
- Data can be transferred between HT monitors or from the HT-10
- Unit can be used as a standalone monitor by itself and data later transferred to a HT monitor
- Standard parameters include - ECG, NiBp, SpO2, Temp, Resp



Pacemaker Intravenous Intracardiac

Pacing monitor to be used with intravenous intracardiac monitoring



HT- Pacemaker Module

- Esophagus pacing
- Non invasive
- Module to be used with HT monitors

Pictures are for illustrative purposes only and may vary from actual product.



Resuscitator Silicone

Infant
Child
Adult



Resuscitator Ambubags

Original Adult
Child Infant



Resuscitator Peep Valves

2 - 10 cm
5 - 20 cm



Resuscitator Peep valve Diverters



Rebreathing Bags



Oxygen Masks



Infant Resuscitator

Easy to use
Gas powered
Designed to supply safe, consistent and optimal resuscitation for infants



Nebuliser Ultrasonic and Piston

Ultrasonic breakdown of molecules
Adjustable vapour settings
Nebuliser rate > 3ml/min
Top of the range



Oxygen Generators

Low noise
Oxygen concentration > 90%
Oxygen Flow 0.5 ~ 3 L/min or 0.5 ~ 5L/min
Output pressure 30 ~ 70 kpa
Options: Flowmeter , monitor Humidifier



Oxygen Cells

Various oxygen cells
For Anaesthesia and Respiratory



Oxygen Cylinders, Bags & Trolleys




Oxygen Monitor Oxiquant

Microprocessor Controlled
Continuous monitoring of Oxygen
Large display
Flexible alarms



Oxygen Blenders

Blender for Precise oxygen
Concentrations between 21 and 100% with a single control knob



Oxygen Generators Plant for Hospitals and Clinics

Molecular Sieve Oxygen generator plant
No more Liquid Oxygen or Bottled Oxygen deliveries
From 150 bed to unlimited size
Backup systems running parallel - if one system down then the other will continue to supply
93 % + Oxygen purity
4 to 5.5 bar output (400 to 550 kpa)

Pictures are for illustrative purposes only and may vary from actual product.



Baby Scales Digital

Digital baby scales with gram
Indication
Up to 20 Kg



Baby Scales Dial

Mechanical baby scales with kg/
pound indication
Up to 20 kg




Table-top Microscale

Auto zero & auto tare
Easy to read display
Precision readability up to 1\15000
Keyboard with audible feedback
Lightweight and durable design
Overload protection for load-cell
Audible beep alerts user when
reaching preset target weight
Plate size: 240 x 290mm
Rechargeable Battery / Mains Oper-
ated



Balance

NHB Series
The NHB – accurate, fast and versa-
tile % weighing function.
Specifications:
ABS Plastic housing, Stainless Steel
top Plate, LCD display with backlight
Underweighing (NHB600 /
NHB6000)
Auto Zero tracking, Audible alarm
for pre-set tare weight, Accumulation
function ' Windshield



Scale VET

Stainless steel platform with rubber
mat
Rechargeable battery / Mains
operated
Platform size: 900 x 500mm
Capacity: Large :150kg x 20g
Small : 60 kg x 10g
Optional: Indicator wall mount
bracket




Adult Scale Mechanical

Die-cast beam
Easy to read from either side
Makes weighing simple and fast
Height rod: 75-192cm
Platform size: 370 x 770mm
Capacity: 180kg x 100g
Castors for easy movement
Brackets or trolleys



Digital Physician Scale

Fully electronic digital scale
Bright 20mm LED display
Auto zero after weight removed
Rechargeable battery / Mains operat-
ed Battery life: 20 hours constant use
Auto shut off to save battery life
Hold function to maintain mass read-
ing Height rod: 810 - 2100mm
Platform size: 275 x 375mm
Capacity: 300kg x 100g



Digital Physician Scale with BMI

BMI (Body Mass Index) Reading
Fully electronic digital scale
Auto zero after weight removed
Rechargeable battery / Mains operat-
ed .Battery life: 20 hours constant use
Auto shut off to save battery life
Hold function to maintain mass read-
ing .Height rod: 810 - 2100mm
Platform size: 275 x 375mm
Capacity: 300kg x 100g



Measuring tape

2 meter height measuring tape
Easy installation



Bathroom scales

Large LCD Display with backlight
Max. Capacity 180kg x 100g
Weighing Modes: kg/lb/st
Platform Size: 320 x 310 x 8mm
Power: 2 x CR2032 lithium batteries

Pictures are for illustrative purposes only and may vary from actual product.



Portable Suction

New technology Oil free lubricating pump
AC / DC 220V Ac and 12V DC
Suction : 15L/Min
Reservoir : 1000ml
Low noise



Suction Unit Mobile

220V 50 Hz low noise Oil-free
2500 ml twin bottles
Low noise
Casters
Suction > 20L / min
Foot switch included



Suction Unit Mobile

220V 50 Hz low noise Oil-free
2500 ml twin bottles
Low noise
Casters
Suction > 20L / min
Foot switch included



Obstetric Suction Cup Extractor

Smooth, soft, malleable silicone ensures optimal adhesion
Smooth external surface with small longitudinal ridges



Silicone Tubing Medical Grade

ID 6.4mm OD 11.1mm
ID 8mm OD 13mm

High grade consistent wall thickness
100 m rolls available



Suction Filters

Bacterial filters for suction units
To fit silicone tubing ID 6.4mm



Wall Suction

Units with regulator
Optional with liners



Stomach Washer

Power Voltage: AC220V 50Hz
Washing Cycle: <40s
Input power: ≤250VA
Washing capacity: 250~350mL/time
Suction capacity: 350~450mL/time
Flux: ≥2.0 L/min
Pressure control: flush and suction pressure set into 47~67kPa, the continue running tolerance no more than ±5kPa



Suction Bottles

Various suction bottles
Glass and PVC

Disposable bottles with liners available



Suction Parts


Piston pumps
Control boards
Battery backup pc boards
Tubing
Shockabsorbers

Pictures are for illustrative purposes only and may vary from actual product.



MS400 Multiparameter Simulator

12 - Lead ECG
Respiration
Temperature
2 - Channel IBP
35 Arrhythmias
3.5" TFT touch screen
Rechargeable Li-Ion battery



MS200 NiBp Simulator

Built in air chamber for Adult and Neonate
Pressure gauge testing, Leak test and Relief valve testing
NiBp simulations for both Arm and wrist
Adult and neonate with eight patient variations



MS100 SpO2 Simulator

Embedded R - Curves
Separate connection between simulator and host
TFT 320 X 240 screen
Rechargeable Li-ion battery



BC401 Urine Analyzer

Testing of GLU, BIL, SG, KET, BLD, PRO, URO, NIT, LEU, VC, PH, MAL, CR and UCA.
High Luminance white LED
2.4" LCD display
Compatible with 8, 10, 11, 12, 14 Strips
USB and bluetooth



Electrosurgical Tester - ESpro

mV, mV Peak, mA, Crest Factor & Wattage displays
Internal data storage for 3 full tests
Large graphical display and cursor selection
Digital data output via USB / RS232



Electrosurgical Tester - ESprime

Capture, Store, Print RF waveforms
Industry Standard RF Current Measurement **for multiple pulsed waveforms**
Automated User-defined testing sequences
USB / Rs232 and Ethernet comm ports
VGA colour touch screen



Defibrillator Analyser DFbase or DFpro

Quick selection of ECG signal types
Defibrillator information retrieval
Monophasic and Biphasic testing
Test for shock algorithms
Measurement of outgoing energy
Fully AED compatible



Infusion pump Analyser - IF+

Simultaneous operation of two chambers
+- 1% of reading flow accuracy
Digital calibration
Large graphic display
Auto chamber recognition



Multi Parameter NiBp

NiBp, Ecg, Temp, IBP
Synchronised IBP output
Flash programmable
ECG output with full NSR waveform
Digital pressure envelope offset
+- 1% of reading pressure accuracy
Small, portable, lightweight



Electrical safety tester SB

Testing for electrical safety of operating equipment
Predefined, standardised test steps
Colour touch screen
Mains independent operation
Internal data memory
Data backup with flash card
USB data interfaces with PC etc.



UV light tester

For testing of Phototherapy
UV light intensity and wavelengths



Bilirubin tester

Neonatal Transcutaneous Jaundice meter

Is used in the dynamic clinical examination of neonate jaundice, by determining the concentration of bilirubin transcutaneously.

Pictures are for illustrative purposes only and may vary from actual product.



Ultrasound Tablet
Chison Sonotouch 30

A dedicated tool designed for vascular access operators
Small footprint, compact and easy to move
LED screen, sharp image
Waterproof touch screen
Dual transducers



Ultrasound S9

15" LCD 45 degree adjustable angle
Smart full touch panel
Two transducer sockets
Linear, convex, micro-convex, Endocavity, Phased, Intraoperative, TEE, Bi-plane, Pencil, Volumetric, Endocavity 4D and Laproscope probes available
Advanced application technology
Built in battery for 90 minutes scanning



Ultrasound 4D colour
HY8000PRO

19" high definition LCD monitor with wide viewing angle
10" touch screen with 15 deg adjustable angle
Height and position adjustable control panel with muted adjustment
5 transducer sockets and an additional socket for CW probe
Additional Endocavity probe holder and gel warmer
Linear, convex, micro-convex, Endocavity, Phased, Intraoperative, TEE, Bi-plane, Pencil, Volumetric, Endocavity 4D and Laproscope probes available
Advanced application technology : TDI, Stress Echo and Elastography
Full patient database solutions: Dicom 3.0, AVI/JPG, USB2.0, HDD, DVD, PDF report



Chison ECO - 1
B & W

- Ultra light—6 kg
- Advanced Image technology: THI, SRA, i-Image, Compound imaging
- Trapezoidal mode



Chison ECO - 3
Expert

- 30 deg rotatable LED screen
- Advanced Image technology: THI, SRA, i-Image, Compound imaging
- Great Image quality : Support B,2B,4B,B/M, M, PW
- Chroma



Chison Q5
Compact colour doppler

- 15 " LCD screen
- Advanced 4D technology: 4T (Fast, Light, Quiet, Smart)
- Professional OB report package
- Trapezoidal Mode
- Advanced Technology : SRA, Compound imaging, THI, i-Image



Chison i3
Compact colour doppler

- 19 " LCD screen, 4 Transducer Con.
- Optional Advanced 4D technology: 4T (Fast, Light, Quiet, Smart)
- Professional OB report package
- Optional Dicom 3.0
- Advanced Technology : SRA, Compound imaging, THI, i-Image



Printers

Sony and Mitsubishi Video printers



Accessories

- Paper : SONY / MITSUBISHI
- Gel in 5 Liters or 250 ml

Pictures are for illustrative purposes only and may vary from actual product.



**Vascular Doppler
BV-520**

Detection of blood flow Arterial / venous at 8Mhz
 Blood Flow velocity
 LED display



**Vascular Doppler
BV-520T**

Unidirectional
 Detection of blood flow Arterial / venous
 Blood flow velocity
 LCD with Waveform



**Vascular Doppler
BV-520T-B**

Bidirectional
 Realtime bloodflow detection
 Blood flow velocity
 Colour LCD with waveform
 Rechargeable battery



**Vascular Doppler
BV-620V**

Bidirectional
 Detection of bloodflow Arterial / venous
 Blood flow velocity
 LCD with Waveform
 Rechargeable Battery



**Vascular Doppler
BV-620VP**

Bidirectional
 Detection of bloodflow Arterial / venous
 Blood flow velocity
 LCD with Waveform
 Rechargeable Battery
 Built in printer – waveform printer

Dopplers - Fetal



**Fetal Doppler
Sonotrax - A Sonotrax - B (display)**

Sonotrax handheld Doppler can meet routine examination requirements of obstetricians, midwives, and pregnant women. It has manual heart rate processing mode. It has audio output, and can connect with earphone or recorder with audio input. Sonotrax has a slim, power saving and interchangeable ultrasound probe.
 A = No display B = Display (Pulse , battery)



**Fetal Doppler
Sonoline C**

45*825mm LCD display
 FHR Performance
 FHR Range: 50-240BPM.
 Waveform display
 Probe Frequency: 2.0MHz
 Working Frequency: 2.0MHz 10%
 Resolution: 1BPM
 Accuracy: 2BPM
 Auto Shut-OFF: After 1 minute no signal, power off automatically.



Pictures are for illustrative purposes only and may vary from actual product.



Batteries

We supply original and repacked batteries.



BP Aneroids

Blood pressure aneroids with bulb, valves and gauge



BP Aneroids Wall

Complete with cuff, bulb and tubing

BP Aneroid Mobile

Complete mobile aneroids



BP Arm Automatic

Automatic measurement of Systolic, Diastolic and Pulse
80 memory
Auto power OFF



BP Wrist Automatic

Automatic measurement of Systolic, Diastolic and Pulse
Auto power OFF



First aid Kits

Regulation 3 / 7
Home kits
Vehicle / motorist kit
Basic life support
Intermediate / Advanced and Nurses kits



Furniture

Beds
Cupboards
Bins
Lockers
Trolleys



HB and Glucometers

Hemoglobin meters
Glucometers
Strips


Mattress

Anti Decubitus
size:2000*980mm
Pressure range:15 - 20KPa



Mattress Foam

Supply all sizes
Incubator sizes available



Needle Destroyer

Syringe needle destroyer



Paper

Recording paper for CTG's
ECG recorders
Ultrasonds
Vital signs monitors
ETC.



Peakflow Meters

Adults
Mouthpieces for lungfunction units



Plaster Saws

12000 revs/min
Speed setting 1-6
Suitable to cut plastic
Blades available

Patient Transfers

For easy and safe movement of patients between beds / trolleys



Resuscitators

Neonatal
Child
Adult
Silicone and rubber
Peep valves/ masks
etc



Vein finders

Akicare SureView

Compact & portable
Clear image
Accurate 3D positioning



Stethoscopes

From nurses to high quality Cardiology stethoscopes



Temperature Infrared

Measurement in 1 sec.
NO touch - forehead
0.1 degree accuracy
Easy to use



Ultrasonic Cleaners

Easy cleaning of surgical instruments, glasses, jewelry
Concentrated cleaner available

Pictures are for illustrative purposes only and may vary from actual product.

AED - Disposable defibrillator pads - Various
 Airways , Guedel
 Alcohol swabs
 Antiseptic solutions
 Bandage clips
 Bandage adhesive elastic
 Bandage conforming
 Bandage crepe
 Bandage cohesive
 Bandage WOW
 Blood pressure cuffs - neonatal
 Body bags with or without windows
 Bowie Dick test packs
 Burn dressings
 Catheters
 Cervical collars
 CPR mouth pieces
 CPR masks
 Eye shields
 First aid kits
 First aid dressings
 Forceps - Budget stainless steel
 Forceps - Budget plastic

Gloves - examination
 Gauze swabs
 Hernia dressings and gauze
 Needle holders
 Plaster roll elastic
 Plaster roll Zinc Oxide
 Plaster strips
 Scalpel handles
 Scalpels
 Scissors budget
 Scissors operate
 Scissors Metzbaum
 SpO2 probes - adult and neonatal
 Syringes
 Wound closure strips
 Vaginal speculums
 Thermometers
 Tongue depressors
 Tracheal dilators
 Umbilical cord clamp
 Urine bags

Service & Repairs



EVOTEC is a CEASA registered corporate member and as such can offer you the highest quality of service .

Evotec follow the ECRI standard of for best service practice and quality repairs

EVOTEC is an Importer and Distributor of medical Test equipment
 We repair and service the following

AED defibs
 Aspirators
 Autoclaves
 Centrifuges
 Defibrillators
 Diagnostic units
 Dopplers
 ECG Recorders
 Electrosurgery units
 Fetal Monitors
 Fetal transducers - Toco & US
 Flow meters

Incubators
 Infant warmers
 Injection pumps
 Monitors
 Nebulisers
 NiBp monitors
 Power supplies
 Phototherapy Lights
 Pulse oximeters
 Radiant warmers
 Regulators
 Scales

Sphygmomanometers
 Spirometers
 SpO2 probes
 Sterilizers
 Suction Units
 Syringe drivers
 Test Equipment
 Temperature monitors
 Tourniquet units
 Ultrasound
 Vital signs monitors
 Warming blanket controllers

Pictures are for illustrative purposes only and may vary from actual product.

EVOTEC

Your one-shop solution
Since 1997

**DIRECT IMPORTER
AND DISTRIBUTOR**

Pictures are for illustrative purposes only and may vary from actual product.

658 Alouette street
Elarduspark
Pretoria

P.O. Box 4333
Rietvalleirand
0174

Phone: 012 345 1176 / 4738 / 6421
Fax : 012 345 4381 or 086 550 4126
E-mail: sales@evotec.co.za

www.evotec.co.za